

Phos-tag application data

Phosphorylation analysis of extracellular signal-regulated kinase (Erk)

~ Comparison of mobility between Bis-Tris-HCl buffer system and Tris-AcOH buffer system (2)

~

● SAMPLE INFORMATION

		MW (kDa)
Protein	Erk1, Erk2 (A431)	200以下
Protein status	normal	-

● ELECTROPHORESIS CONDITION

Gel	8.0% polyacrylamide (Bis-TrisHCl, Tris-AcOH)
Phos-tag conc.	25 μ M
Metal complex	Zn ²⁺

Visualization	immunoblotting
Antibody	anti-pT202/Y204 anti-Erk1/2

● ASSAY FLOW

- 1 EGF stimulation of cell lysate
- 2 Phos-tag electrophoresis
- 3 Immunoblotting

● RESULT

- Phosphorylation status and phosphorylated bands of activated isoform obvious.
- In the two buffer system, phosphorylation sites of separated isoform were identified by anti-Thr-202 antibody and anti-Tyr-204 antibody.

● NOTE

Another data of simultaneous analysis : 007-1 007-3

● REFERENCE

Phos-tag SDS-PAGE systems for phosphorylation profiling of proteins with a wide range of molecular masses under neutral pH conditions. Kinoshita E, Kinoshita-Kikuta E, Koike T. : *Proteomics*, **12**, 192 (2012)

key words : Bis-Tris-HCl, Tris-AcOH, Zn²⁺-Phos-tag, Erk, Erk1, Erk2, A431, HMW